


# Specialised Oral Care Products From **Mobile Medical Systems Intl**

T: 02 9252 4410 / E: info@mobilemedical.com.au

## Surround<sup>®</sup> Toothbrush



**Designed for 'Special Needs', elderly and orthodontic patients**

- 3 headed Surround<sup>®</sup> Toothbrush 'hugs and surrounds' teeth with 45 degree angled bristles to provide simultaneous brushing of gums and the front, back and biting surfaces optimising plaque removal.

- Using simple back and forth brushing motion, more tooth surfaces brushed with each stroke; oral care completed in shorter period of time with teeth cleaned more completely and efficiently over traditional, conventional toothbrushes.

- Outstanding tool for those who need assistance in oral care, have difficulty holding or using a toothbrush.

### Features

- \* Compact head with soft, end-rounded bristles for comfortable and safer brushing
- \* Side bristles positioned at 45 degree angle recommended by the ADA to clean both teeth and gums
- \* *FLEXIBLE* structure and soft tactile surface makes gripping easier and helps protect oral environment from trauma in sudden unexpected movement during brushing.
- \* Covered in a non latex, rubbery material with padded head and neck.
- \* Head & neck of brush designed to make brushing safer even with those who tend to bite down on the toothbrush



**Now in Adult and Toddler Size**

