

Rainbow® Stabilising Systems*

Size guides: board and wrap



Mobile Medical
Systems

RSS minimises sudden, unsafe patient movement during procedures; complies with **AAPD Guidelines on Protective Stabilisation for Paediatric Dental Patients** that devices be “*easily used, appropriately sized for patient, soft & contoured to minimise potential injury, specifically designed for patient stabilisation (i.e. not improvised equipment) and able to be disinfected.*”



Rainbow® Boards

Height	Width	Rainbow® Board	Dark Teal	Cream
2 Feet	11"	Infant Board	RBLdt	RBLc
3 Feet	11"	Small Board	RBSdt	RBSc
4 Feet	12"	Medium Hinged Board	RBHMDt	RBHMc
4 Feet	12"	Medium Straight Board	RBMdt	RBMc
5 Feet	16"	Large Hinged Board	RBHLdt	RBHLc
5 Feet	16"	Large Straight Board	RBLdt	RBLc

What Size to Buy?

Up to this approx. weight/age:

Recommendation

26 pounds.....	(Infants).....	Infant Board - Infant Wrap
35 pounds.....	(1-4 years old).....	Small Board - Small Wrap
60 pounds.....	(5-8 years old).....	Medium Board - Medium Wrap
145 pounds.....	(9-12 years old).....	Large Board - Large Wrap
185 pounds.....	(Teen/ Adults).....	Large Board - X-Large Wrap
250 pounds.....	(Teen/ Adults).....	Large Board - XX-Large Wrap

Recommended Board and Wrap sizes, based on approximate age and weight of patient.

Infant Wrap 50.8 cm length / 109.22 cm wide

11.8 kg / Infants

Small Wrap 63.5 cm length / 88.9 cm wide

15.9 kg / 1 – 4 years old

Medium Wrap 81.3 cm Long / 111.8 cm wide

27.2 kg / 5 – 8 years old



*MIPS- Medical Immobilisation and Protective Stabilisation